# OPHTHAMIC GAS / RETINAL GAS

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ARCEOLE	SF6	C2F6	C3F8
Retinal Detachment (pure gas for retinopexy)	✓		
Diabetic Retinopathy	✓		
Retinal Detachment (mix gas/air after a vitrectomy)		✓	
Macular Hole		✓	
Giant Tears	✓	✓	
PVR			✓
Ocular Trauma	✓	✓	
			

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<b>ARCEOLE</b>	 SF6	 C2F6	 C3F8
Expansion delay	1 day	1.5 day	3 days
Effective tamponade time (half of the eye volume)	6 days	12 days	30 days
Presence of gas	15 days	30 days	60 days
Rate of expansion	2	3.3	4
Concentration for a non-expansive gas	20 %	16-17 %	12 %
Maximum volume of gas	1.5 ml	1 ml	0.8 ml
			

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


Poorly soluble in water

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### ADVANTAGES

Guarantees the presence in the eye for a ± long term.

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Chemical inertness

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**No toxic**, low interaction with organic tissue.

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Easy visualization of the posterior segment.

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## 특징 CHARACTERISTICS

High difference of surface tension  
with water (71 mN/m at 20°C)

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Very low density/water  
낮은 밀도/물



## 장점 ADVANTAGES

Effective sealing of tears  
효과적 열공 폐쇄

Creation of a unique bubble  
→ *Filling whole space*  
→ *Good cover, even if the retina is folded*

High reapplication force on the tissue  
조직에 재적용력


→ *progressive reattachment of the retina*  
→ *drainage of sub-retinal fluids*  
→ *stabilisation of the retina*

**PRECAUTION** : avoid a massive displacement of sub-retinal fluid in the macula through manipulation of the globe

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


Initial expansion induces an increased pressure. Progressive dissolution induces a decreasing tamponade surface.

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

\*The favourite gas injection for simple RD : 0.3ml of SF<sub>6</sub>

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Sub-retinal passage

- Haemorrhage
- Choroidal /sub-conjunctival passage

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Hypertony

- Retinal tear
- Cataract
- Inflammation
- Endophthalmitis
- Keratopathy

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A patient implanted with intraocular gas must absolutely

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- Avoid pressure variations

(elevation, travel by plane, diving with or without isobar caisson, hyperbar caisson treatment)

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A patient who is **not able to respect head position** (obese, pregnant)

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1. Lockwood AJ, Yang YF. Nitrous oxide inhalation anaesthesia in the presence of intraocular gas can cause irreversible blindness. Br Dent J. 2008 Mar8;204(5):247-8

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- Easy to inject
- No secondary surgery
- 3 gas for short-medium term tamponade
- Strong reapplication force
- Efficient sealing
- Disappear gradually

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- Simple RD
- Macular Hole,
- Diabetic retinopathy
- Traumatic Retinal Detachment

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- Stay in the eye
- Permanent tamponade
- Long term sealing
- Preservation of the surface tamponade
- Difficult to inject and remove
- No reapplication pressure
- Potential emulsification

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- RD with proliferation
- Giant Retinal Tears ( $\geq 180^\circ$ )
- RD with AIDS

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Atmospheric pressure  
prefilled syringe

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Perfluoropropane (C<sub>3</sub>F<sub>8</sub>)

**LIQUIFIED GAS UNDER PRESSURE**

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

NOTE: Do not use the contents after the expiration date printed on the cylinder label.

**DESCRIPTION**

ISPAN\* Perfluoropropane (C<sub>3</sub>F<sub>8</sub>) is a liquified gas under pressure and is administered by injection into the vitreous cavity. It is Octafluoropropane (C<sub>8</sub>F<sub>18</sub>) from the Haloalkanes chemical family. The boiling point is -36.7°C (-34.1°F) and the vapor pressure at 20°C is 100 psig (pounds per square inch gauge). Perfluoropropane is clear and colorless with a faintly sweet odor. ISPAN\* C<sub>3</sub>F<sub>8</sub> purity: perfluoropropane (Octafluoropropane) 99.9% (minimum), air 1000 ppm (maximum), and perfluoropropene 10 ppm (maximum).

**INDICATIONS**

ISPAN\* Perfluoropropane (C<sub>3</sub>F<sub>8</sub>) is a surgical aid for use in the treatment of uncomplicated retinal detachment by pneumatic retinopathy. It is used in the form of an intravitreal injection for selected retinal breaks and to aid in resorption of subretinal fluid. Associated measures used include transconjunctival and transcleral cryotherapy and laser photocoagulation.

**CONTRAINDICATIONS**

Proliferative vitreoretinopathy (PVR) greater than Stage C, the mental or physical inability to maintain the therapeutic position for 5 postoperative days, severe glaucoma with more than a minimum of field loss and a cup to disc ratio equal to or greater than 0.6; uveitis; severe peripheral retinal degeneration; and high altitude travel, included but not limited to airline travel.

**MODE OF ACTION**

During the healing phase, the surface tension of the gas can prevent further progression of the retinal detachment by holding the retina against the choroid and permitting the retinal pigment epithelial pump to remove the subretinal fluid. The perfluoropropane diffuses from the eye in approximately 6 to 8 weeks.

**WARNINGS**

Use of Nitrous Oxide (N<sub>2</sub>O) must be stopped at least 10 minutes before gas injection to ensure an adequate postoperative bubble is achieved. Do not administer Nitrous Oxide (N<sub>2</sub>O) if a gas bubble is present. Nitrous Oxide (N<sub>2</sub>O) rapidly partitions into the gas bubble causing expansion and a pressure increase in the eye that has been known to result in vision decrease and blindness. There is a risk of cataract formation if the lens is inadvertently damaged by the needle during gas injection during pneumatic retinopathy.

Do not use the ISPAN\* gas if the cylinder pressure is below 50 psi as the expansion performance of the gas may change resulting in elevated intraocular pressure.

Acute rises in intraocular pressure (IOP) which threaten ocular blood flow for greater than 10 minutes should be controlled with paracentesis of aqueous fluid or removal of part of the gas bubble. Patients with compromised ocular blood flow such as those with severe diabetic retinopathy or ocular ischemia are at greater risk of vascular occlusion following the use of an expansile gas bubble. The intraocular pressure (IOP) should be checked by an experienced surgeon with either tactile touch or applanation tonometry when ISPAN\* C<sub>3</sub>F<sub>8</sub> is in place. Schiøtz tonometry will give false low values compared to the true IOP.

Patient positioning following intravitreal gas injection is of great importance. The bubble must be properly situated with the proper positioning to allow contact of the gas bubble against the retinal hole or holes internally. Prone or seated face down positioning can prevent protracted contact between the gas bubble and the lens to avert a posterior subcapsular cataract, as well as to prevent pressure on the ciliary body and iris, and to prevent pupillary block in aphakic patients, which might increase intraocular pressure. The central retinal artery should be monitored during and after gas injection. Administration of systemic carbonic anhydrase inhibitors or topical glaucoma medications may be given for less severe elevations of intraocular pressure.

- Air travel is contraindicated until the gas bubble has completely dissipated. Normal cabin pressure changes will cause a severe enlargement of the gas bubble with a resultant increase in IOP.<sup>1,2</sup>
- Patients should not travel through high elevations and over mountain ranges until the gas bubble has dissipated.<sup>3</sup>
- Patients should not receive hyperbaric oxygen therapy until the gas bubble

**Do not use the ISPAN\* gas if the cylinder pressure is below 50 psi as the expansion performance of the gas may change resulting in elevated intraocular pressure.**

are provided with this product and should be given to the patient prior to discharge from their eye surgery. The patient card is a convenient way to remind the patient about the important restrictions noted above, including limitations on the use of Nitrous Oxide (N<sub>2</sub>O) in subsequent surgical or dental procedures, travel in an airplane or through high elevations; and when to remove the bracelet. The bracelet is to be worn by the patient to alert subsequent health professionals that the patient may have a gas bubble in their eye and to confer with the ophthalmologist prior to treating the patient. Ensure both sides of the bracelet and card are completed and reviewed with the patient. Additional cards and bracelets may be obtained from Alcon Customer Service at 1-800-862-5266, or your local Alcon representative.  
Reorder: Hospital bracelets - TAM110; Patient Cards - TAM111

1. Lincoff, H, Weinberger, D, Reppucci, V, Lincoff, A: Air travel with intraocular gas. I. The mechanisms for compensation. Arch Ophthalmol 107:902-906, 1989
2. Lincoff, H, Weinberger, D, and Stergiu, P: Air travel with intraocular gas. II. Clinical considerations. Arch Ophthalmol 107:907-910, 1989
3. Hanscom, TA, and Diddle, KR: Mountain travel and intraocular gas bubbles, AM J Ophthalmol 104:546, 1987
4. Jackman, SV, and Thompson, JT: Effects of hyperbaric exposure on eyes with intraocular gas bubbles. Retina 15:160-166, 1995

**PRECAUTIONS**

Caution should be used in eyes with angle recession, pigment dispersion syndrome, significant anterior synechiae, traumatized eyes and eyes with significant vitreous hemorrhage obscuring an adequate view of the peripheral retina.

Sterile surgical techniques should be used for injection of ISPAN\* C<sub>3</sub>F<sub>8</sub>. Blepharitis or other lid infections should be treated as if for intraocular surgery prior to using ISPAN\* C<sub>3</sub>F<sub>8</sub> for pneumatic retinopathy. Endophthalmitis has been reported rarely following pneumatic retinopathy.

There are no known teratogenic effects of ISPAN\* C<sub>3</sub>F<sub>8</sub> when injected into the eye. Until such information is available, it should be used with caution in pregnant women.

Sterility cannot be assured when the gas is transferred from the tank to a sterile syringe. The gas must be filtered through a sterile 0.22 µm filter prior to injection into the eye and used immediately. A pressure reducing gas regulator should be used to remove ISPAN\* C<sub>3</sub>F<sub>8</sub> from the cylinder. The delivery pressure of the gas should not exceed 10 psig. The lecture bottle stand is recommended for maintaining the necessary upright position of the gas cylinder during use. Close cylinder valve when not in use.

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

Plastic syringe is not a right solution because of air permeation through a plastic syringe-container at atmospheric pressure.

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
Because of air permeation, manufacturer increases the initial concentration of gas by 50% over concentration mentioned in the IFU

**Ready to use!**

**EasyGas SF<sub>6</sub>**  
G-80950  
20% Sulfurhexafluoride  
with 80% air,  
Syringe 40 ml, sterile

EasyGas SF <sub>6</sub>	
Effective tamponade time [days]	6
Retention time / Longevity* [weeks]	1-2
Non-expansive gas concentration* [%]	20

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- **Arceole<sup>®</sup> by Arcadophta** (SF<sub>6</sub>; C<sub>2</sub>F<sub>6</sub>; C<sub>3</sub>F<sub>8</sub>) patented<sup>®</sup> 특허보유  
30mL single dose
- **Arceole<sup>®</sup> Set Multidose by Arcadophta** (SF<sub>6</sub>; C<sub>2</sub>F<sub>6</sub>; C<sub>3</sub>F<sub>8</sub>)  
30mL single dose
- **MicroGas<sup>®</sup> by Micromed** (SF<sub>6</sub>; C<sub>2</sub>F<sub>6</sub>; C<sub>3</sub>F<sub>8</sub>) patent infringement 특허침해  
50mL single dose
- **GasMate<sup>®</sup> by Howard Instruments** (SF<sub>6</sub>; C<sub>2</sub>F<sub>6</sub>; C<sub>3</sub>F<sub>8</sub>) patent infringement 특허침해  
50mL single dose
- **GOT<sup>®</sup> by Alchimia** (SF<sub>6</sub>; C<sub>2</sub>F<sub>6</sub>; C<sub>3</sub>F<sub>8</sub>) patent infringement 특허침해  
75mL single dose

ARCEOLE Pressurized gas

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VS.

Liquefied gas

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## ARCEOLE

### Low-pressurized gas (가압가스)

1 to 5 uses => **short duration storage**

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Small containers => **easy storage**

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1 set sterile => **easy handle by the surgeon only**

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1 set with accessories => **no additional cost, manipulation by the surgeon only**

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Pressurized container => **No toxic residue**

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### Liquefied gas (액화가스)

1 to 500 uses => **risk of waste** (shelf life, tightness)

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**manipulation by two nurses + 1 surgeon**

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Need of sterile complementary components (tube, filter, syringe, cannula, gas distributor) & investment on non sterile intermediate equipment (manometer, Regulator, feet..) + biomedical maintenance =>

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**Raise of toxic residues** with gasconsumption Multiples uses

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		Low-pressurized gas 가압가스	Liquefied gas 액화가스
Additional safety 추가적인 안전성	Usage out surgical theatre 수술장 밖의 사용법	Yes (Retinopexy)	No
	Complete traceability of the device per patient 환자별 추적성	Yes Complete set delivered with traceability stickers	No (bottle – multi use of the same can) Syringe traceability?
	Ease of the follow-up of the shelf-life of the consumable 용이한 유통기한 확인	Yes Shelf-life on the package	Hard No real control
Legal responsibility 법적 책임	Transport 운송	No constraint 1000 units by air	Restricted transport forbidden by air
	Legal cover of the marking EC Responsibility of the Industrial 마킹 EC의 법적 표지 산업의 책임	Everything including can, gas, mixture device, sterilization process and use in ophthalmic surgery	Tightness of the bottle Composition of gas
	Responsibility of the surgeon / hospital 외과 의사 / 병원의 책임	Surgery only	Selection of the medical device for preparation and surgery

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		Low-pressurized gas 가압가스	Liquefied gas 액화가스
Economics 경제적	Initial investments 초기 투자	Very low (stock)	Important Double pressure gauges & reducer of pressure, tripod
	Purchase of consumables 소모품구매	No (complete set)	Yes (tube, syringe, filter, needle)
	Maintenance 유지비	No	Yes (each year)
	Optimal availability of gas stock according to indications 적응증에 따른 최적의 가스 재고 확보	Low investment SF6/C2F6/C3F8	Huge investment (1 to 3 available gas)
	Time mobilization of staff for preparation of consumables 소모품 준비를위한 직원의 시간 동원	Low	Important
	Return on investment (Surgeon, private hospital, hospital) 투자수익	High	Low
	Economic loss 경제적 손실	Low / Rare (when shelf-life have been passed)	High volume of non consumed gas (flushing ...)

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